



Marketing **100+** Countries / Regions

*PUSEN Endoscope  
R&D and Production Base*

**PUSEN**

Zhuhai Pusen Medical Technology Co., Ltd.  
5/F, Building 1, No 33, Keji San Road,  
High-tech Zone, Tangjiawan Town,  
519085 Zhuhai, Guangdong, People's Republic of China  
(86) 400-688-2023

Rev.202401

**PUSEN**

**Single Use  
Cysto-Nephro Videoscope  
PC200**





**Integrated Suction Port**

Clear vision, rapid drainage

Custom Function Buttons

**Innovative Patented Handle**

Ergonomic Handle to relieve fatigue

**Integrated Straight Working Channel**

Straight-Line instrument entry, making the operation more concise

**Kink-Resistant Tube**



**Dual LED**

Clear vision for accurate diagnosis and treatment

**HD Imaging**

1280\*720 HD image algorithm, more precise during treatment

**In-Package Connecting Tube**

Connecting the irrigation port with the irrigation system

**Larger Deflection Angle**

No blind spots in endoscopic examination, enhancing diagnostic and therapeutic effectiveness

**Fully-Flexible Lens Body**

Alleviating patient discomfort

**Suction Control**



Smooth cystoscopy without frequent catheter changes  
Controlled bladder filling reduces patient discomfort  
Shorten cystoscopy time, improve diagnosis of early and small lesions

15Fr

**Specification**

Field of View	90°
Diameter of Insertion Portion	15 Fr
Minimum Instrument Channel Width	7.2Fr
Bending Angle	U/D 210°